

# Squidly Session

GUIDE




## User Guide

# Joining a Squidly Meeting

Click the secure meeting ID or copy it into your preferred internet browser.

- 1 Make sure camera, microphone and sound options are toggled on.
- 2 If the connection doesn't form after 10 seconds try refreshing as the host/participant.

Sessions are designed for one-to-one connections only.



## Why can't I join a Squidly meeting?

- 1 Ensure you have a stable internet connection.
- 2 Contact the host to make sure the meeting ID is valid.
- 3 Try on a different browser. Google chrome is advised.
- 4 Contact Squidly support at [contact@squidly.com.au](mailto:contact@squidly.com.au).




# Using Squidly Webcam Eye-Gaze

Tap access in the menu and select calibrate.


- 1 To re/calibrate tap calibrate.
- 2 Follow the on-screen prompts on the calibration window and carefully follow the guide with your eyes.

Once a calibration is complete you will receive an accuracy score.



If the eye-gaze cursor doesn't appear immediately after calibration it may be turned off.

- 3 Tap access in the menu and toggle the eye-gaze on/off; or
- 4 Go to settings in the menu and make sure host/participant eye-gaze is enabled/disabled.




# Eye-Gaze Calibration

Tap settings in the menu and navigate to host/participant calibration.


1 Increase/decrease calibration grid size.

2 Increase/decrease calibration speed.

3 Change the calibration guide image.




## Improve Eye-Gaze Accuracy



## Using the AAC Grid

Tap access in the menu and select aac.


- 1 Interact together using the aac default board, 'voco chat'; or
- 2 Search for an aac board.
- 3 Tap speak to say the full sentence and clear to delete the sentence.



## Starting a Quiz or Assessment


Tap share in the menu and select quiz.

- 1 Select or search for a quiz and begin the quiz together.
- 2 Navigate through quiz questions and answers and tap submit when you have completed the quiz.
- 3 Receive a quiz score and download an AI-summary report once you've submitted.




# Settings

Squidly enables you to change both host and participant settings. Tap settings in the menu.



## Accessible Games and Tools

Tap share in the menu and select apps. Choose from our growing catalogue of accessible games and tools designed to keep clients motivated and engaged in therapy and learning sessions.



# Access Features

Squidly is designed to work seamlessly with a wide range of access methods, ensuring that every user can engage comfortably and independently. The platform automatically adapts to your personal setup and device configuration.

- 1 Integrated eye-gaze access - works with any compatible eye-gaze system set to Windows Control.
- 2 Switch scanning - supports multiple switch inputs (keyboard: space/backspace).
- 3 Touch access - fully touch-responsive interface for tablets or touch-enabled devices.
- 4 Browser-based - runs directly in modern web browsers with no software installation required.
- 5 Text-to-speech - converts written text into clear, natural, custom speech outputs.
- 6 Keyboard and mouse compatibility - works with standard or adaptive input devices for flexible access.
- 7 Adjustable interfaces - customisable visual layouts, colour contrast, languages, and pointer dwell times.
- 8 Host navigation control - allows the host to navigate the entire session on behalf of the participant.

# Meeting Controls

